β-ketonen teststrips + code-chip ketonen

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TaiDoc Technology Corporation
B1-7F, No.127, Wugong 2nd Rd.,
Wugu Dist., 24888 New Taipei City, Taiwan

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For use with TD-4289B

Use

Thanks to the Diatesse XPER β-ketone test strips, you can use the Diatesse XPER meter (TD-4289B) to measure your β-ketone levels or have them checked by healthcare professionals. The system uses both fresh capillary whole blood samples extracted from a finger and venous whole blood. For *in vitro* diagnostic use only (not for internal use). The system is not designed to diagnose and screen for diabetes mellitus. Healthcare professionals can use the test strips to test both capillary and venous whole blood; only capillary whole blood samples can be tested at home.

Code chip (calibration) instructions

Calibration

Recalibrate the meter every time you open a new box of β-ketone test strips by inserting the code chip in your meter.

Contact points

Calibration instructions


- Insert the code chip (the section with the contact points, with the text facing up) into your deactivated meter. Wait until the code number and "KET" appear on the display.
- Hoove the ketone code chip. "OFF" will appear on the display and the meter will switch itself off.

Test results cannot be accurately displayed if the code on your meter does not match the calibration code on the box of β-ketone test strips.

Checking the code

Check whether the correct code (it must match the code on the box of β-ketone test strips you are using) and "KET" are displayed on the meter's display before performing a test. If the codes do not match, do not perform the test, but instead input the correct code chip (again). If the codes still do not match, then please contact our customer service department.

Description of the β-ketone test strips

- Absorbent opening**
This opening absorbs a drop of blood. As soon as this part of the β-ketone test strip comes into contact with your blood, it automatically absorbs it.
- Control window**
The window allows you to check whether the absorbent opening of the β-ketone test strip has absorbed sufficient blood.
- Grip**
Hold this section when inserting the β-ketone test strip into your meter's slot.
- Contact points**
Insert this section into your meter's slot. Ensure that the β-ketone test strip is firmly inserted.

Your β-ketone tests

Wash and dry your hands before conducting a test.

Step 1
Insert the β-ketone test strip (the section containing the contact points) firmly into your meter's slot. When the β-ketone test strip is correctly inserted, the meter will perform a number of self-tests.

Step 2
Take a blood sample using the β-ketone test strip. In order to provide an accurate test result, sufficient blood must be drawn. Allow the absorbent opening of the β-ketone test strip to touch the drop of blood and wait until the control window is filled. Do NOT apply the strip to a sample that has been smeared. Once the control window contains sufficient blood, the meter starts an automatic countdown.

Step 3
After a few seconds, the meter will display your β-ketone levels. The last measurement is automatically stored in the meter. Switch the meter off by removing the test strip. Dispose of the used β-ketone test strip.

For further information, please refer to the user manual for the meter.

The used lancet and the β-ketone test strip could constitute a potential biological hazard. Carefully dispose of them in accordance with the local regulations.

Reading your results

Your β-ketone measurements are plasma equivalent values and are shown in millimoles of β-ketones per litre of blood (mmol/L) or in milligrams of β-ketone per deciliter of blood (mg/dL).

The β-ketone test measures beta-hydroxybutyrate (BHB), the most important of the three β-ketones in your blood. Under normal circumstances, the concentration of β-ketones in your blood should be less than 0.6 mmol/L (6.2 mg/dL).¹¹

The β-ketone concentration can increase if you are fasting and/or dieting, exercising intensively or suffer from diabetes and fall ill. If your β-ketone test results in a "Lo", repeat the test using a new test strip. If you receive the same result or if the result does not match how you are feeling, then contact your healthcare professional. Always follow the recommendations of your healthcare professional before making any changes to your diabetes treatment. If the levels for your β-ketone test are between 0.6 and 1.5 mmol/L (6.2 and 15.6 mg/dL), this could be due to an issue requiring professional medical attention. Carefully follow your healthcare professional's advice. If your β-ketone test levels are greater than 1.5 mmol/L (15.6 mg/dL), immediately contact your healthcare professional, as this means you could develop diabetic ketoacidosis (DKA).

¹¹: Wiggam M,I., O’Kane M.J., Harper R., Atkinson A.B., Hadden D.R., Trimble E.R., Bell P.M. Treatment of diabetic ketoacidosis using normalization of blood 3-hydroxybutyrate concentration as the endpoint of emergency management. Diabetes Care 1997; 20: 1347-52.

Consult your doctor in order to determine the target range that suits you best.

Dubious or inconsistent results

- If your test results are unusual or do not match the way you are feeling:
- Ensure that the β-ketone test strip's control window is filled with blood.
 - Check the expiry date of the β-ketone test strips.
 - Check that your meter and the β-ketone test strip are working properly using the control solution.

PLEASE NOTE: Unusually high or low β-ketone values could be symptoms of a serious medical condition. If the majority of your test results are unusually high or low, contact your healthcare professional.

Chemical components

β-Hydroxybutyrate dehydrogenase (*Pseudomonas sp.*) ≥ 0,5 U

Mediator 55%

NAD ≥ 0,5 µg

Enzyme-protector 8%

Non-reactive ingredients 29%

Quality control tests

The control solution contains a fixed quantity of β-ketones that react with the β-ketone test strips. You can check your meter, the β-ketone test strip and the testing technique you use by comparing the control solution to the range printed on the label of the β-ketone test strip packaging. Please see the meter's user manual for complete instructions for testing the device.

▲ The reference range of the control solution can vary between bottles of control solution and between packages of β-ketone test strips, so always check the range printed on the label of the β-ketone test strip packaging you are currently using.

Supplementary information for healthcare professionals

Always wear gloves and follow your institution's infection prevention procedures when conducting tests involving blood samples. Only use fresh blood samples. Healthcare professionals can use the β-ketone test strips to test capillary and venous whole blood samples.

Sample size: 0.8 µl

Reaction time: 10 seconds

Measurement range: 0.1 - 8.0 mmol/L (1.0 - 83.2 mg/dL)

HCT range: 10 - 70%

Warnings

- For *in vitro* diagnostic use (not for internal use).
- Single use only.
- Healthcare professionals and other users that test multiple patients using this system must treat all objects that have been in contact with human blood with utmost care, including sterilised objects, in order to prevent the transmission of infectious diseases.
- Read this package insert and the user manual for your Diatesse XPER meter before using these test strips. Only use **Diatesse XPER β-ketone test strips** with the Diatesse XPER meter (TD-4289B) in order to obtain accurate results and for the device to remain under the manufacturer's warranty.
- Results could be inaccurate if tests are conducted on patients with abnormally low blood pressure or in a state of shock.
- For patients with poor peripheral circulation, we do not recommend taking a capillary blood sample from the approved sampling locations, as the results may not accurately reflect the physiological β-ketone levels. This could be the case in the following circumstances: serious dehydration due to diabetic ketoacidosis or due to stress-induced hyperglycaemia, hyperglycaemic hyperosmolar non-ketotic coma, shock, NYHA Class IV heart failure or peripheral arterial occlusive disease (PAOD).
- Test strips and lancets must always be kept out of reach of children. If swallowed, consult a doctor immediately.

Restrictions

- Hematocrit:* The HCT level is restricted to values between 10 and 70%.
- In vitro* ascorbic acid up to 228 µmol/L, captopril up to 23 µmol/L, cholesterol up to 15 mmol/L, dopamine up to 5.8 µmol/L, genticis acid up to 117 µmol/L, L-DOPA up to 10 mg/L, paracetamol up to 1324 µmol/L, triglycerides up to 30 mmol/L, uric acid up to 3 mmol/L and unconjugated bilirubin up to 400 µmol/L will not significantly influence the results of a test.
- Altitude:* An altitude up to no more than 3,500 meters shall not influence the test results.

Storage and use

- ▲ Do not use the test strips after their expiry date**
 - Store the β-ketone test strips in a cool and dry place at a temperature between 2 °C and 30 °C (35.6 °F and 86 °F) with a relative humidity of between 10% and 85%.
 - Keep the β-ketone test strips away from direct sunlight.
 - Always keep the β-ketone test strips in their original foil packaging. Do not unwrap a β-ketone test strip from its foil packaging until it is required.
 - Do not touch the β-ketone test strips with wet hands.
 - Use the β-ketone test strip immediately after removing it from the foil packaging.
 - Except when removing a β-ketone test strip, always keep the foil packaging closed.
 - Do not bend, cut or alter the β-ketone test strip.

Accuracy

The reference method is β-Hydroxybutyrate LiquiColor®. The reagents can quantitatively measure the presence of β-ketones that leads to the suspicion of a patient suffering from diabetic ketoacidosis.

β-ketones	Capillary samples	Range, mean
Regression	$y = 0.922x + 0.1863$ ($y = 0.9193x + 1.7787$)	

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