Bloedglucose-teststrips

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TaiDoc Technology Corporation
B1-7F, No.127, Wugong 2nd Rd.,
Wugu Dist., 24888 New Taipei City, Taiwan

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For use with

Warnings

- For *in vitro* diagnostic use only (not for internal use).
- Single use only.
- Healthcare professionals and other users that test multiple patients using this system must treat all objects that have been in contact with human blood with utmost care, including sterilised objects, in order to prevent the transmission of infectious diseases.
- Read this package insert and the user manual for your Diatesse XPER meter before using these test strips. Only use **Diatesse XPER blood sugar test strips** with the Diatesse XPER meter (TD-4289B) in order to obtain accurate results and for the device to remain under the manufacturer's warranty.
- Results could be inaccurate if tests are conducted on patients with abnormally low blood pressure or in a state of shock.
- For patients with poor peripheral circulation, we do not recommend taking a capillary blood sample from the approved sampling locations, as the results may not accurately reflect the physiological blood sugar levels. This could be the case in the following circumstances: serious dehydration due to diabetic ketoacidosis or due to stress-induced hyperglycaemia, hyperglycaemic hyperosmolar non-ketotic coma, shock, NYHA Class IV heart failure or peripheral arterial occlusive disease (PAOD).
- Test strips and lancets must always be kept out of reach of children. If swallowed, consult a doctor immediately.

Use

Thanks to the Diatesse XPER blood sugar test strips, you can use the Diatesse XPER meter (TD-4289B) to measure your blood sugar levels or have them checked by healthcare professionals. The system uses fresh capillary whole blood samples extracted from a finger and venous whole blood. The system is not designed to diagnose and screen for diabetes mellitus. Healthcare professionals can use the test strips to test both capillary and venous whole blood; only capillary whole blood samples can be tested at home.

Restrictions

- Hematocrit*: The HCT level is restricted to values between 0 and 70 L/L.
- Metabolites*: Dopamine, L-Dopa, methylropa, uric acid, ascorbic acid and paracetamol at a normal blood concentration shall not significantly influence the results of a blood sugar test.
- No significant interference has been observed in blood sugar tests when there is a presence of galactose, maltose or fructose.
- Lipaeamic effects*: Triglycerides in the blood up to 33. 9 mmol/L(3,000 mg/d L) do not significantly influence the results, but can have an impact at higher levels.
- Altitude*: Altitude up to no more than 3,275 meters shall not influence the test results.
- Only use heparin as an anticoagulant for fresh capillary or venous whole blood.
- There are **no** restrictions/interferences with respect to maltose up to 1,000 mg/dL and icodextrin up to 2,000 mg/dL.

The following substances can, when they are in excess of the value limits, return increased blood sugar values when testing using the Diatesse XPER meter (TD-4289B):

<i>Name</i>	<i>Limiting concentration (mg/dL)</i>	<i>Therapeutic/physiological concentration range (or upper limit) mg/dL</i>
Acetaminophen (paracetamol)	> 20	0.45 - 3
Ascorbic acid	> 5	2
Uric acid	> 10	2 - 8

Storage and use

PLEASE NOTE: Do not use the blood sugar test strips after their expiry date

- Note the date that you first opened the tube of blood sugar test strips on the label (only for strips in a tube). Do not use the blood sugar test strips after their expiry date.*
- Store the blood sugar test strips in a cool and dry place at a temperature between 2 °C and 30 °C (35.6 °F and 86 °F) with a relative humidity of between 10% and 90%.
- Keep the blood sugar test strips away from direct sunlight. Do not keep the blood sugar test strips in very humid locations.
- Always keep the blood sugar test strips in their original tube/foil packaging. Do not place the blood sugar test strips in a different tube. When the blood sugar test strips are contained in foil packaging, do not unwrap them until required.
- Use the blood sugar test strip immediately after removing it from the tube/foil packaging. Reseal the tube immediately after removing a blood sugar test strip.
- Except when removing a blood sugar test strip, always keep the tube/foil packaging closed.
- Do not bend, cut or alter the blood sugar test strip.

Description of the blood sugar test strips

- Absorbent opening**
Apply a drop of blood at this point. As soon as this part of the blood sugar test strip comes into contact with your blood, it automatically absorbs it.
- Control window**
The window allows you to check whether the absorbent opening of the blood sugar test strip has absorbed sufficient blood.
- Grip**
Hold this section when inserting the blood sugar test strip into your meter's slot.
- Contact points**
Insert this section into your meter's slot. Ensure that the blood sugar test strip is firmly inserted.

PLEASE NOTE

The front of the blood sugar test strip must be facing upwards when you insert it into the meter. If the contact points are not wholly and firmly inserted into the meter's slot, your test results could be incorrect.

Your blood sugar tests

Wash and dry your hands before starting.

- Step 1**
Insert the blood sugar test strip (the section containing the contact points) firmly into your meter's slot. When the blood sugar test strip is correctly inserted, the meter will perform a number of self-tests.
- Step 2**
Take a blood sample of 0.5 µl using the blood sugar test strip. In order to provide an accurate test result, sufficient blood must be drawn. Allow the absorbent opening of the blood sugar test strip to touch the drop of blood and wait until the control window is filled. Do NOT apply the strip to a sample that has been smeared. Once the control window contains sufficient blood, the meter starts an automatic countdown.
- Step 3**
After a few seconds, the meter will display your blood sugar levels. The last measurement is automatically stored in the meter. Switch the meter off by removing the test strip. Dispose of the used blood sugar test strip.

For further information, please refer to the user manual for the meter.

The used lancet and the blood sugar test strip could constitute a potential biological hazard. Carefully dispose of them in accordance with the local regulations.

Reading your results

Your blood sugar measurements are *plasma equivalent values* and are shown in millimoles of glucose per litre of blood (mmol/L). The measuring range of the meter is from 0.55 to 44.44 mmol/L (10-800 mg/dL).

<i>Time of day</i>	<i>Normal plasma glucose range for nondiabetics</i>
Before meals on an empty stomach	< 5.6 mmol/L (100 mg/dL)
2 hours after meals	< 7.8 mmol/L (140 mg/dL)

Source: American Diabetes Association. Standards of medical care in diabetes - 2018 Jan; 41 (Supplement 1): S1-S2.

Consult your doctor in order to determine the target range that suits you best.

Dubious or inconsistent results

If your test results are unusual or do not match the way you are feeling:

- Ensure that the blood sugar test strip's control window is filled with blood.
- Check the expiry date of the blood sugar test strips.
- Check that your meter and the blood sugar test strip are working properly using the control solution.

PLEASE NOTE: Unusually high or low blood sugar values could be symptoms of a serious medical condition. If the majority of your test results are unusually high or low, contact your healthcare professional.

Quality control tests

The control solution contains a fixed quantity of glucose that reacts with the blood sugar test strips. You can check your meter, the blood sugar test strip and the testing technique you use by comparing the control solution to the range printed on the label of the blood sugar test strip tube/foil packaging. Please see the meter's user manual for complete instructions for testing the device.

PLEASE NOTE: The reference range of the control solution can vary between bottles of control solution and between packages of blood sugar test strips, so always check the range printed on the label of the blood sugar test strip tube/foil packaging you are currently using.

Chemical components

Glucose dehydrogenase (E. coli)	8%
Electron-shuttle	55%
Enzyme-protector	8%
Non-reactive ingredients	29%

Supplementary information for healthcare professionals

Always wear gloves and follow your institution's infection prevention procedures when conducting tests involving blood samples. Only use fresh blood samples. Healthcare professionals can use the blood sugar test strips to test capillary and venous whole blood samples.

Sample size: 0.5 µl

Reaction time: 5 seconds

Measurement range: 0.6 - 44.44 mmol/L (10 mg/dL- 800 mg/dL)

HCT range: 0 - 70 L/L

Accuracy

Diabetes experts have suggested that blood sugar meters must be within ± 0.83 mmol/L (± 15 mg/dL) of the reference method when the glucose concentration is less than 5.55 mmol/L(100 mg/dL) and within ± 15% of the reference method when the glucose concentration is 5.55 mmol/L(100 mg/dL) or higher.

<i>Reagents system LOT</i>	<i>Within ±0.28 mmol/L</i>	<i>Within ± 0.56 mmol/L</i>	<i>Within ± 0.83 mmol/L</i>
TD16J11S-CEE	44 / 58 (75.9%)	56 / 58 (96.6%)	58 / 58 (100%)
TD16K122-CEE	45 / 58 (77.6%)	58 / 58 (100%)	58 / 58 (100%)
TD16L112-CEE	40 / 58 (69.0%)	58 / 58 (100%)	58 / 58 (100%)

For a glucose concentration of < 5.6 mmol/L, 100% of the measured results were within ± 0.83 mmol/L.

On average, 98.9% were within ± 0.6 mmol/L.

<i>Reagents system LOT</i>	<i>Within ± 5%</i>	<i>Within ± 10%</i>	<i>Within ± 15%</i>
TD16J11S-CEE	100 / 142 (70.4%)	133 / 142 (93.7%)	140 / 142 (98.6%)
TD16K122-CEE	102 / 142 (71.8%)	138 / 142 (97.2%)	142 / 142 (100%)
TD16L112-CEE	100 / 142 (70.4%)	139 / 142 (97.9%)	142 / 142 (100%)

For a glucose concentration of ≥ 5.6 mmol/L, an average of 96.2% of the measured results were within ± 10%.

Note: When the results of the Diatesse XPER blood sugar test strips are compared to the reference values, the difference values beneath 5.55 mmol/L(100 mg/dL) are expressed in mmol/L, while the difference values above 5.55 mmol/L(100 mg/dL) are expressed in percentages.

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